



# FULL LINE PRODUCT CATALOG

## BE THE DIFFERENCE®

PDI offers a broad range of evidence-based, market-leading *Interventional Care*, *Environment of Care*, and *Patient Care* solutions, all designed to help reduce preventable infections, control healthcare costs, and, ultimately, help save lives. It's all part of our passion to *Be The Difference*, every day.

## Environment of Care

The environment is a primary, and often unexpected, source of infections. Maintain a clean and sanitary environment with our *Environment of Care* solutions, the #1 disinfecting wipes in healthcare.\*

### Sani-Cloth® Prime Germicidal Disposable Wipes and Sani-Prime® Germicidal Spray



Powered by ACCELOQUAT™ formulation, protect your patients, staff, and facility from HAIs in just 1 minute.



### Super Sani-Cloth® Germicidal Disposable Wipe



Ideal for daily use in fast-paced environments that require a short contact time and broad coverage of microorganisms.



### Sani-Cloth® AF3 Germicidal Disposable Wipe



Alcohol free. Fragrance free. Worry free. Ideal for use around patients and staff, especially those with respiratory sensitivities.



### Sani-Cloth® Plus Germicidal Disposable Cloth



A low-level disinfecting wipe for general surface disinfection.



### Sani-Cloth® Bleach Germicidal Disposable Wipe



Ideal for disinfecting high-risk areas endemic with *Clostridium difficile* spores and norovirus.



### Easy Screen® Cleaning Wipe

A compatible touchscreen cleaner that works at the speed of technology.



- Features the power of 70% isopropyl alcohol (IPA)
- Fast-drying, no residue, anti-fogging, anti-streaking
- Rapidly cleans dirt, grime, fingerprints, and smudges

\*GHX 2017 Q4 trend report

Different infection risks require different solutions. Our *Interventional Care* solutions help you protect even the most vulnerable patients.

## Prevantics® Skin Antiseptics

3.15% (w/v) Chlorhexidine Gluconate (CHG) and 70% (v/v) IPA topical skin antiseptic.

- Swabstick and Maxi Swabstick are approved for preoperative skin preparation and are effective on dry site with as little as 30 seconds of scrubbing and 30 seconds of dry time
- Swabstick and Maxi Swabstick are clinically proven to have continuous antimicrobial activity for up to 7 days
- Swab approved for preinjection skin preparation
- Premoistened, ready-to-use applicators contain no glass



## Prevantics® Device Swab

3.15% (w/v) CHG/70% (v/v) IPA solution for the disinfection of needleless access sites prior to use.

- Confirmed efficacy at the 5-second scrub time<sup>1</sup> and 5-second dry time
- Strip format hangs conveniently on an IV pole for point-of-care accessibility



## Profend™ Nasal Decolonization Kit

Easy-to-use, pre-saturated 10% PVP-iodine swabstick – just snap and swab!

- Kills 99.7% of *S. aureus* at 1 hour and 99.9% at 12 hours<sup>2</sup>
- Ideal for *S. aureus*- and MRSA-colonized patients in the surgical unit, ICU, and other areas of the hospital



**Reference: 1.** Hayden MK, et al. A Randomized Cross-Over Clinical Trial to Compare 3.15% Chlorhexidine/70% Isopropyl Alcohol (CHG) vs. 70% Isopropyl Alcohol Alone (Alcohol) and 5s vs. 15s Scrub for Routine Disinfection of Needleless Connectors (NCs) on Central Venous Catheters (CVCs) in an Adult Medical Intensive Care Unit (ICU). Oral Abstract Presented at 2014 ID Week Conference, October 11, 2014, Philadelphia, PA. **2.** PDI *in vivo* Study 0113-CTEVO.

# Patient Care

Help your patients and residents remain clean and comfortable with our *Patient Care* solutions.

## Sani-Hands® Instant Hand Sanitizing Wipes

Designed for everyone – including staff, patients, residents, and visitors.

- Contains 70% Ethyl Alcohol, Aloe & Vitamin E emollients, and is clinically proven to kill 99.99% of germs<sup>1</sup>
- FDA food code compliant. Safe to use before eating and drinking when used as directed<sup>2</sup>
- Proven to remove dirt and soil better than gels<sup>3,4</sup>



## Hygea® Multi-Purpose Washcloths

- Suitable to use on all areas of the body, including face and neck
- Ideal for perineum use and incontinence care
- Gentle cleansing formula contains soothing Aloe & Vitamin E, hypoallergenic, alcohol-free and latex-free



## Hygea® Flushable Personal Cleansing Cloths

- Premoistened washcloth for full body personal cleansing and incontinence care
- Cleans and refreshes better than dry bathroom tissue alone
- Hypoallergenic, alcohol-free and latex-free



**References:** 1. Clinical data is representative of Sani-Hands® (NDC #: 10819-50) 65.9% Ethyl Alcohol formulation. 2. Food Code Recommendations of the United States Public Health Service Food and Drug Administration, 2013. 3. Clinical data is representative of Sani-Hands (NDC #: 10819-50) 65.9% Ethyl Alcohol formulation vs. the Purell (NDC #: 21749 515) 62% Ethyl Alcohol formulation. 4. Data on file, PDI, Orangeburg, New York.

**Products kill bugs, but people prevent infections.**™ Our point-of-care accessories help put PDI products in the right place at the right time to help healthcare professionals do their job more efficiently and effectively.

All of PDI's point-of-care accessories are designed to support compliance with facility protocols.

## 3-in-1 Universal Sani-Bracket®

Easily attach **Sani-Cloth** and **Sani-Hands** canisters to walls or equipment for convenient access. Customize with product-specific labels and adaptors to meet your facility's particular needs.



## Sani-Canister Caddy®

The countertop holder for canisters is perfect for nursing stations and carts. It helps display important product and compliance information, such as contact time and precautionary language.



## Sani-TAG® Equipment ID System

**Sani-TAG** Equipment ID System adheres to surfaces and equipment to help staff identify the specific surface disinfectant that should be used on the particular item.

The "Clean Tags" help staff easily identify when equipment has been disinfected, and meet regulatory and accreditation standards for identifying clean vs. dirty equipment.

All **Sani-TAG** labels remove easily and leave no sticky residue.



## PDI Safety Stand

The stand holds one canister of product and has optional attachments that can be affixed directly to the floorstand to hold personal protective equipment.



# Reorder Information

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI

## ENVIRONMENT OF CARE



### Sani-Cloth® Prime Germicidal Disposable Wipe & Sani-Prime™ Germicidal Spray

Powered by ACCELOQUAT™ formulation, protect your patients, staff, and facility from HAIs in just 1 minute.

Large Canister	P25372	6 x 6.75"	12/160s	16.00 lbs	1.410 ft	10/3
X-Large Canister	P24284	7.5 x 15"	6/70s	27.84 lbs	1.010 ft	10/4
Large Individual Packet	H06182	5 x 8"	10/50s	7.99 lbs	0.540 ft	20/4
X-Large Individual Packet	U13195	11.5 x 11.75"	3/50s	8.65 lbs	0.500 ft	24/3
Spray Bottle	X12309	n/a	9/32oz per bottle	19.24 lbs	0.828 ft	12/4



### Super Sani-Cloth® Germicidal Disposable Wipe

Fast, 2-minute contact time allows for a quick room turnover. Ideal for daily use in fast-paced environments that require short contact times and broad coverage of microorganisms.

Large Canister	Q55172	6 x 6.75"	12/160s	26.00 lbs	1.55 ft	10/3
X-Large Canister	Q86984	7.5 x 15"	6/65s	15.30 lbs	1.01 ft	10/4
Large Individual Packet	H04082	5 x 8"	10/50s	7.72 lbs	0.539 ft	20/4
X-Large Individual Packet	U87295	11.5 x 11.75"	3/50s	8.00 lbs	0.50 ft	24/3



### Sani-Cloth® AF3 Germicidal Disposable Wipe

Protecting your staff, patients, and residents has never been easier! Alcohol and fragrance-free formula is ideal for use around those with respiratory sensitivities.

Large Canister	P13872	6 x 6.75"	12/160s	26.00 lbs	1.546 ft	10/3
X-Large Canister	P63884	7.5 x 15"	6/65s	15.30 lbs	1.01 ft	10/4
Large Individual Packet	H59200	5 x 8"	10/50s	7.72 lbs	0.54 ft	20/4
X-Large Individual Packet	U27500	11.5 x 11.75"	3/50s	8.00 lbs	0.50 ft	24/3
Portable Pack	M928S80	8.2 x 9.8"	9/100s	22.63 lbs	0.95 ft	10/4
Pail	P1450P	7.5 x 15"	2/160s	12.40 lbs	0.95 ft	10/5
Refill (for Pail)	P2450P	7.5 x 15"	2/160s	11.10 lbs	0.53 ft	18/4

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI



### Sani-Cloth® Plus Germicidal Disposable Cloth

Ideal for use in alternate care settings, including physician and dental offices.

Large Canister	Q89072	6 x 6.75"	12/160s	26.00 lbs	1.546 ft	10/3
X-Large Canister	Q85084	7.5 x 15"	6/65s	15.30 lbs	1.012 ft	10/4



### Sani-Cloth® Bleach Germicidal Disposable Wipe

Ideal for disinfecting high-risk areas endemic with *Clostridium difficile* spores and norovirus.

Large Canister	P54072	6 x 10.5"	12/75s	26.25 lbs	1.54 ft	10/3
X-Large Canister	P25784	7.5 x 15"	6/65s	20.50 lbs	1.01 ft	10/4
Large Individual Packet	H58195	5 x 7"	10/40s	8.10 lbs	0.51 ft	19/5
X-Large Individual Packet	U26595	11.5 x 11.75"	3/40s	8.20 lbs	0.50 ft	24/3
Pail	P7007P	7.5 x 15"	2/160s	16.50 lbs	0.95 ft	10/5
Refill (for Pail)	P700RF	7.5 x 15"	2/160s	15.20 lbs	0.53 ft	18/3



### Sani-Cloth® HB Germicidal Disposable Wipe

Designed to accommodate the many situations and hard, nonporous environmental surfaces found in healthcare settings.

Large Canister	Q08472	6 x 6.75"	12/160s	26.00 lbs	1.55 ft	10/3
X-Large Canister	Q85484	7.5 x 15"	6/65s	15.30 lbs	1.01 ft	10/4



### Easy Screen® Cleaning Wipe

A compatible touchscreen cleaner that works at the speed of technology. Contains 70% IPA.

Medium Canister	P03672	6" X 9"	12/70s	10.64 lbs	.990 ft	8/5
-----------------	--------	---------	--------	-----------	---------	-----

## INTERVENTIONAL CARE



### Profend™ Nasal Decolonization Kit

Uniquely designed swabsticks pre-saturated with 10% Povidone-Iodine (PVP-I) — just snap and swab.

Patient Kit	X12048	n/a	4 swabs/ patient pack, 12 patient packs/shelf unit, 4 shelf units/case	2.70 lbs	0.263 ft	35/5
-------------	--------	-----	---	----------	----------	------

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI



**Prevantics® Skin Antiseptics**

3.15% (w/v) CHG and 70% (v/v) IPA skin antiseptic.

Swab	B10800	n/a	10/100s	4.30 lbs	0.27 ft	21/7
Swabstick	S40750	n/a	10/50s	7.30 lbs	0.98 ft	15/3
Maxi Swabstick	S41950	n/a	10/30s	12.90 lbs	1.53 ft	5/5
Swab (for kits)	B11400	n/a	1/3000s	18.93 lbs	1.02 ft	10/4
Swabstick (for kits)	S32450	n/a	1/500s	5.80 lbs	1.02 ft	10/4
Compact Swabstick (for kits)	S42850	n/a	1/500s	6.60 lbs	0.65 ft	10/8
Maxi Swabstick (for kits)	S27350	n/a	1/300s	13.20 lbs	1.02 ft	10/4



**Prevantics® Device Swab**

For disinfecting needleless access sites prior to use. Contains 3.15% (w/v) CHG and 70% (v/v) IPA solution.

Swab	B19600	n/a	10/100s	4.30 lbs	0.28 ft	21/7
Swab Strip	B123ST	n/a	10/160s	4.68 lbs	0.52 ft	8/5



**Alcohol Prep Pads - Sterile**

Provides rapid-acting skin antiseptics in a sterile formulation. 70% IPA.

Medium Packet	B60307	1.1 x 2.6"	20/200s	12.40 lbs	0.84 ft	8/7
Large Packet	C69900	1.7 x 3.5"	10/100s	7.20 lbs	0.48 ft	12/8



**PVP Iodine Prep Pads and Swabsticks**

Saturated with a 10% USP Povidone-Iodine prep solution. A broad-spectrum antiseptic for patient pre-operative skin preparation.

Medium Packet	B40600	1.1 x 2.6"	10/100s	4.61 lbs	0.22 ft	24/8
Large Packet	C12400	1.7 x 3.5"	10/100s	10.23 lbs	0.52 ft	12/7
Swabstick – 1s	S41350	n/a	10/50s	7.50 lbs	0.87 ft	13/4
Swabstick – 3s	S41125	n/a	10/25s	8.70 lbs	0.87 ft	13/4



**PVP Iodine Scrub Swabsticks**

Saturated with a 7.5% Povidone-Iodine cleansing solution. A broad-spectrum antiseptic for scrubbing skin prior to surgery.

Swabstick – 1s	S48050	n/a	10/50s	7.50 lbs	0.87 ft	13/4
Swabstick – 3s	S82125	n/a	10/25s	8.30 lbs	0.87 ft	12/4



**PVP Iodine Duo-Swab® Prep & Scrub Swabsticks**

Contains one cleansing scrub swabstick and one antiseptic swabstick in a convenient 2-step package.

Swabstick – 1s	S23125	n/a	10/25s	7.13 lbs	0.87 ft	13/4
----------------	--------	-----	--------	----------	---------	------



REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI

## PATIENT CARE

**Sani-Hands® Instant Hand Sanitizing Wipes**

Contains 70% Ethyl Alcohol, Aloe &amp; Vitamin E emollients, and is clinically proven to kill 99.99% of germs.

Medium Canister	P13472	6 x 7.5"	12/135s	18.45 lbs	0.99 ft	8/5
Large Canister	P15984	6 x 7.5"	6/220s	15.25 lbs	0.80 ft	10/5
Individual Packet	D43600	5 x 8"	10/100s	9.86 lbs	0.65 ft	9/7
Bedside Pack	P71520	5.5 x 8.4"	48/20s	11.20 lbs	0.46 ft	15/6

**Hygea® Multi-Purpose Washcloths**

Strong, disposable washcloths that are ideal for cleansing and refreshing patients.

Premium - Solo® Softpak	J14143	11.5 x 9.5"	6/60s	11.00 lbs	0.67 ft	9/7
Solo Softpak	J14108	12 x 8"	9/60s	13.30 lbs	0.67 ft	10/5
Solo Softpak	J21396	12 x 8"	6/96s	13.30 lbs	0.72 ft	10/6
Resealable Softpak	J22750	12 x 8"	12/48s	13.30 lbs	0.67 ft	10/6
Bulk Individual Packet	U12095	11.5 x 10"	1/400s	16.70 lbs	1.56 ft	6/4

**Hygea® Flushable Personal Cleansing Cloths**The ideal flushable pre-moistened washcloth for personal general skin cleansing and incontinence care. Safe for use with sewer and well-maintained septic systems.<sup>a,b,c</sup>

Solo Softpak	A500F48	5.3 x 6.8"	12/48s	6.89 lbs	0.37 ft	15/7
--------------	---------	------------	--------	----------	---------	------

**Hygea® BZK Antiseptic Towelette (alcohol free)**

For use wherever antiseptic cleansing/refreshing is needed. Contains 0.13% Benzalkonium Chloride.

Individual Packet	D35185	7 x 5"	20/100s	17.30 lbs	1.26 ft	9/3
-------------------	--------	--------	---------	-----------	---------	-----

**Hygea® Obstetrical Towelette**

For personal hygiene at every napkin change and whenever cleansing/refreshing is needed. Saturated with Benzalkonium Chloride, Alcohol 5%, Chlorothymol.

Individual Packet	D74800	7.75 x 5"	10/100s	9.67 lbs	0.65 ft	9/7
-------------------	--------	-----------	---------	----------	---------	-----

<sup>a</sup>When used as directed<sup>b</sup>Based on the testing guidelines to determine flushability set by the Association of the Nonwoven Fabrics Industry (INDA)<sup>c</sup>Not for use with basement pump systems. Only flush one wipe at a time

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI



**Hygea® Sterile Saline Wipes**

Helps reduce the risk of cross contamination. Convenient, cost-effective alternative to bottled saline and gauze pads. Delivers a consistent dose of sterile saline with every application.

Individual Packet	C22370	6 x 4"	24/24s	13.37 lbs	0.49 ft	23/4
-------------------	--------	--------	--------	-----------	---------	------



**Castile Soap Towelette**

For general skin cleansing, catheter maintenance and clean catch urine samples.

Individual Packet	D41900	7.7 x 5"	10/100s	12.00 lbs	0.65 ft	9/7
-------------------	--------	----------	---------	-----------	---------	-----



**Nice 'n Clean® Baby Wipes**

Developed for babies' soft and delicate skin and offering a variety of products for every stage of baby care. Available in scented, unscented, and sensitive formulas.

Scented - Solo® Softpak	M225XT	7.9 x 6.6"	12/80s	12.70 lbs	0.75 ft	15/4
Scented - Resealable Softpak	Q34540	7.9 x 6.6"	12/40s	6.39 lbs	0.38 ft	15/8
Unscented - Solo Softpak	M233XT	7.9 x 6.6"	12/80s	12.70 lbs	0.75ft	15/4
Unscented - Resealable Softpak	Q70040	7.9 x 6.6"	12/40s	6.39 lbs	0.38 ft	15/8
Newborn Sensitive - Resealable Softpak	Q387S40NY	7.9 x 6.6"	12/40s	6.39 lbs	0.38 ft	15/8



**Sterile Lubricating Jelly**

Facilitates the entry of medical and therapeutic devices. Intended for medical purposes. Will not affect surgical instruments, rubber or plastics. Water soluble, bacteriostatic, latex-free, non-staining and greaseless.

Individual Packet – 2.7 gm	T00137	n/a	12/144s	13.78 lbs	1.18 ft	10/4
Individual Packet – 5.0 gm	T00250	n/a	12/72s	11.40 lbs	0.90 ft	12/4



**Adhesive Tape Remover Pads**

Removes adhesive tape from skin. Reduces patient discomfort.

Individual Packet	B16400	1.2 x 2.6"	10/100s	3.60 lbs	0.25 ft	20/6
-------------------	--------	------------	---------	----------	---------	------



**Io-Gone® Iodine Remover Pads**

Removes extraneous iodine stains from skin and clothing. Quickly dissolves residual iodine color.

Individual Packet	B47000	1.2 x 2.6"	10/100s	4.30 lbs	0.33 ft	17/8
-------------------	--------	------------	---------	----------	---------	------



**Electrode Skin Prep Pads**

Cleanses skin in preparation for electrode placement. Pumice removes excess skin oil for better electrode-to-skin contact and longer interface-free contact.

Individual Packet	B59800	1.2 x 2.6"	10/100s	4.00 lbs	0.33 ft	17/8
-------------------	--------	------------	---------	----------	---------	------

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI



### Nail Polish Remover Pads

For the removal of nail polish in a professional or hospital use. Does not contain acetone.

Individual Packet	B71200	1.2 x 2.6"	10/100s	5.00 lbs	0.33 ft	17/8
-------------------	--------	------------	---------	----------	---------	------



### See Clear® Lens Cleaning Wipes

Ideal for cleaning eye protection gear and glasses. The anti-fog, anti-static, and anti-streaking formula dries quickly.

Individual Packet	D25341	6 x 5"	12/120s	9.59 lbs	0.89 ft	10/5
-------------------	--------	--------	---------	----------	---------	------

## FOOD SAFETY



### Sani Professional® Cleaning Multi-Surface Wipes

Ideal for use on hard surfaces, such as dining tables and countertops.

Solo® Softpack (75 count)	A972FW	9 x 7"	20/75s	14.42 lbs	0.986 ft	10/3
Solo® Softpack (90 count)	A580FW	7 x 11.5"	12/90s	24.00 lbs	0.87 ft ft	10/3

## POINT-OF-CARE ACCESSORIES



### 3-in-1 Universal Sani-Bracket®

Attach **Sani-Cloth** and **Sani-Hands** canisters to walls or equipment for convenient access; customize with product specific labels.

3-in-1 Universal Sani-Bracket (10 pack)	P010900	n/a	10	19.50 lbs	2.03 ft	4/4
3-in-1 Universal Sani-Bracket (1 pack)	P010901	n/a	1	1.75 lbs	0.17 ft	48/5
Sani-Bracket for X-Large Box	P44210	n/a	1	0.70 lbs	0.07 ft	72/8
3-in-1 Sani-Bracket Pole Accessory Kit	P010600	n/a	1	1.65 lbs	0.05 ft	63/14
Bracket Adapter	P015300	n/a	10	0.23 lbs	n/a	n/a



### Sani-Canister Caddy®

Countertop placeholder for canisters. Displays important product and compliance information.

Sani-Canister Caddy Sani-Hands (10 pack)	P013500	n/a	10	1.50 lbs	0.25 ft	28/6
Sani-Canister Caddy Sani-Hands (5 pack)	P013501	n/a	5	0.80 lbs	0.11 ft	30/10
Sani-Canister Caddy Super Sani-Cloth (10 pack)	P012600	n/a	10	1.50 lbs	0.25 ft	28/6
Sani-Canister Caddy Super Sani-Cloth (5 pack)	P012601	n/a	5	0.80 lbs	0.11 ft	30/10
Sani-Canister Caddy Sani-Cloth AF3 (10 pack)	P015200	n/a	10	1.50 lbs	0.25 ft	28/6
Sani-Canister Caddy Sani-Cloth AF3 (5 pack)	P015201	n/a	5	0.80 lbs	0.11 ft	30/10
Sani-Canister Caddy Sani-Cloth Bleach (10 pack)	P013100	n/a	10	1.50 lbs	0.25 ft	28/6
Sani-Canister Caddy Sani-Cloth Bleach (5 pack)	P013101	n/a	5	0.80 lbs	0.11 ft	30/10
Sani-Canister Caddy Easy Screen	P016101	n/a	5	0.80 lbs	0.11 ft	30/10

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI



**Sani-Cloth® AF3 Portable Pack Bracket**

Attach the portable pack to mobile equipment for ease of use and accessibility.

Sani-Cloth AF3 Portable Pack Mobile Equipment Bracket	P013200	n/a	1	0.10 lbs	n/a	n/a
Sani-Cloth AF3 Portable Pack Bracket	P039100	n/a	10	n/a	n/a	n/a



**Sani-Bracket® Deluxe Multi-Canister Floor Stand**

Place in entryways and other high traffic areas. Holds up to 5 canisters.

Sani-Bracket Deluxe Multi-Canister Floor Stand	P013600	n/a	1	41.00 lbs	8.01 ft	2/2
--	---------	-----	---	-----------	---------	-----



**PDI Safety Stand**

Holds one canister of product and has optional attachments that can be affixed directly to the floor stand to hold personal protective equipment.

PDI Safety Stand	P013700	n/a	1	18.00 lbs	1.81 ft	3/7
PPE Station	P014100	n/a	1	5.50 lbs	0.64 ft	11/7



**Sani-TAG® Equipment ID System**

Adhere to equipment to help staff easily identify it as clean. Removes easily and does not leave a sticky residue. Meets regulatory and accreditation standards for identifying clean vs. dirty equipment.

Sani-TAG-Equipment ID System (Clean Tags)	P012200	n/a	8/100s	5.50 lbs	0.29 ft	12/5
---	---------	-----	--------	----------	---------	------

*Sani-TAG - Disinfect with Sani-Cloth Bleach	*Sani-TAG - Disinfect with Sani-Cloth AF3
*Sani-TAG - Disinfect with Sani-Cloth Plus	*Sani-TAG - Disinfect with Sani-Cloth HB
*Sani-TAG - Disinfect with Super Sani-Cloth	*Sani-TAG - Please use Easy Screen Cleaning Wipe

\*Request from sales representative

For more information, or to place an order, contact your local PDI sales representative.

